

OpenSpectrum — Mobile Wireframe Concepts (V1)

Product Direction

Core Principle

The app should feel like:

“A calm memory extension for overwhelmed caregivers.”

The UX prioritizes:

- ultra-low friction
 - fast capture during stressful moments
 - natural-language journaling
 - AI-assisted structuring
 - multiple entries throughout the day
 - longitudinal pattern detection
 - emotional neutrality and support
-

PRIMARY NAVIGATION

Bottom navigation:

1. Home
2. Timeline
3. Quick Log
4. Insights
5. Settings

Floating global action:

- Large microphone button for instant voice capture
-

SCREEN 1 — HOME / TODAY DASHBOARD

Purpose

Immediate awareness + rapid entry.

The parent may open this screen 10–30 times/day.

Mobile Wireframe

OpenSpectrum

Maya ▼

TODAY — Tuesday, May 19

Overall Day So Far

😊

Calm periods: 3

⚠

Dysregulation events: 2

🌙

Sleep last night: 6h 20m

🗣

TAP TO SPEAK

“Quickly describe what happened...”

QUICK LOG

[Sleep] [Food] [Mood] [Medication]

[Meltdown] [Therapy] [Sensory] [School]

RECENT ACTIVITY

2:45 PM — Difficult school transition

1:10 PM — Ate partial lunch

11:30 AM — OT session

8:15 AM — Magnesium supplement

EARLY INSIGHTS

“Sleep under 7 hrs may correlate with increased afternoon dysregulation.”

Home | Timeline | + | Insights | Settings

Purpose

Target completion time:

- 15–45 seconds

AI extracts:

- transition difficulty
- appetite issue
- calming intervention
- supplement timing

Mobile Wireframe

```

| ← Back                                     Recording... |
|-----|
|                                                                          |
|                                     00:18 elapsed                                     |
|                                     ~ waveform visualization ~                                     |
|-----|
| LIVE TRANSCRIPT                                     |
| "Bad transition after school. Refused lunch.                                     |
| Better after weighted blanket..."                                     |
|-----|

```

AI DETECTED	
[Transition Difficulty]	
After School	
[Low Appetite]	
Lunch	
[Helpful Intervention]	
Weighted Blanket	
[Supplement Timing]	
Magnesium — Late	

[Edit] [Save Entry]	

SCREEN 3 — QUICK TAP LOG

Purpose

Ultra-fast structured logging.

Used for:

- repetitive events
- medication
- quick observations
- high-stress moments
- school pickup/dropoff

Goal:
Under 10 seconds.

Mobile Wireframe

Quick Log	

WHAT ARE YOU LOGGING?	
[😊 Good Mood] [😞 Dysregulated]	
[😴 Sleep] [🍽️ Food]	
[💊 Medication] [🧠 Sensory]	
[🗣️ Communication] [🏫 School]	

[👤 Toileting]	[💡 Therapy]
[★ Success]	[📝 General Note]

OPTIONAL DETAILS	
Severity: ☐ Low ☒ Medium ☐ High	
Notes: _____	

[Cancel]	[Save]

SCREEN 4 — DAILY TIMELINE

Purpose

Chronological “life stream.”

This becomes the primary review surface.

Parents naturally think chronologically.

Design Goals

- highly scannable
 - emotionally calm
 - easy filtering
 - visually lightweight
-

Mobile Wireframe

Timeline Filters ⚙️

Tuesday, May 19
7:10 AM
🌙 Sleep ended
6h 20m total sleep

8:15 AM

	🔊 Magnesium	
	Logged by Dad	

	11:30 AM	
	💖 Occupational Therapy	
	"Responded well to sensory activities."	

	2:45 PM	
	⚠️ Difficult school transition	
	Severity: Medium	

	3:10 PM	
	★ Weighted blanket helped recovery	

	6:40 PM	
	🍽️ Ate partial dinner	

	Home Timeline + Insights Settings	

SCREEN 5 — INSIGHTS DASHBOARD

Purpose

Simple, emotionally meaningful pattern detection.

Avoid:

- overwhelming analytics
- medical dashboards
- complex charts early on

Focus on:

- actionable observations
- gentle language
- trend awareness

Mobile Wireframe

	Insights	

	POSSIBLE PATTERNS	

	🌙 Sleep & Regulation	
	"Days with less than 7 hours of sleep	
	showed increased afternoon dysregulation."	

	🧘 Helpful Interventions	
	"Weighted blanket usage often preceded	
	calmer evenings."	

	🏠 Environmental Triggers	
	"Crowded environments may correlate with	
	increased stress markers."	

	📈 Positive Trend	
	"Transition difficulty appears reduced	
	over the past 3 weeks."	

	[View Weekly Summary]	

SCREEN 6 — WEEKLY SUMMARY

Purpose

Low-cognitive-load caregiver reflection.

Useful for:

- parent awareness
- partner communication
- doctor visits
- therapy review

Mobile Wireframe

	Weekly Summary	

	WEEK OF MAY 13-19	

	😊 Best Days	
	Tuesday, Thursday	

	⚠ Most Difficult Periods	

	After-school transitions	
	🌙 Average Sleep	
	7h 05m	
	☆ Most Helpful Supports	
	Weighted blanket	
	Outdoor walks	
	Protein snack	
	📝 Emerging Notes	
	"Late meals may correlate with irritability."	

	[Export Report]	

SCREEN 7 — CHILD PROFILE

Purpose

Support:

- multiple children
- multiple caregivers
- longitudinal profile context

Mobile Wireframe

	Maya Profile	

	Age: 11	
	Communication: Verbal	
	Sensory Profile: Auditory-sensitive	
	CARE TEAM	
	Mom ✓	
	Dad ✓	
	Grandma ✓	
	OT Therapist (view only)	

	COMMON SUPPORTS	
	• Weighted blanket	

• Quiet room	
• Compression vest	
• Protein snacks	

PRIVACY SETTINGS	
Local-first storage ✓	
Encrypted backup ✓	
Anonymous research sharing OFF	

CORE UX PRINCIPLES

1. Multiple Entries Per Day Are Expected

This is not a “daily diary.”

The system should support:

- many small logs
- interrupted workflows
- fragmented observations
- quick emotional unloading

The timeline becomes cumulative intelligence.

2. Voice Is a Primary Interface

Not an add-on.

The app should assume:

- hands may be occupied
- parents may be stressed
- typing may be difficult in the moment

Voice-first dramatically reduces friction.

3. Imperfect Data Is Acceptable

Parents should never feel pressure to:

- fully complete forms
- categorize perfectly
- remember exact times

Partial data over time is still valuable.

4. Emotional Safety Matters

Avoid clinical/judgmental language.

Prefer:

- “needed support”
- “difficult transition”
- “high sensory load”

Avoid:

- “bad behavior”
 - “non-compliant”
 - “problem child”
-

RECOMMENDED MVP PRIORITIES

Build Order

Phase 1

- Home dashboard
- Voice journaling
- Quick log
- Timeline

Phase 2

- AI extraction refinement
- Weekly summaries
- Basic insights

- Export reports

Phase 3

- Caregiver collaboration
- Smart reminders
- Context-aware suggestions

Phase 4

- Anonymous cohort intelligence
- Similar-pattern matching
- Community learning loops

STRATEGIC DIFFERENTIATOR

Most competitors optimize for:

- clinical tracking
- compliance logging
- therapy documentation

OpenSpectrum should optimize for:

- reducing caregiver mental load
- capturing lived experience naturally
- extracting longitudinal meaning
- surfacing useful patterns gently

That distinction is likely the strongest product advantage.